

# Disposable ENT Suction Coagulator

The MEDSTAR's disposable suction coagulators provide precise coagulation with controlled fluid evacuation. The ergonomic handle combined with a flexible, kink-resistant cannula maximizes suction while allowing the surgeon to customize the instrument to the specific procedure at hand.

	○	○	○
French Gauge	12	10	8
Diameter (mm)	4	3.3	2.7
Diameter (inches)	0.158	0.131	0.105

**Specialties**  
Tonsillectomy  
Adenoidectomy

Six inch flexible Shaft with 2-layer insulation may reduce the potential for oral burns



Hand Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord, Single Use

Suction Control : Regulated by placing finger over suction vent on handle



3-prong plug

PE01-8Fr	Hand Control 8Fr	20PCS/Box
PE01-10Fr	Hand Control 10Fr	20PCS/Box
PE01-12Fr	Hand Control 12Fr	20PCS/Box

Foot Control, 6" (15.2 cm) cannula, 10ft. (3.0 m) cord, Single use

Suction Control : Regulated by placing finger over suction vent on handle



Banana plug

PE02-8Fr	Foot Control 8Fr	20PCS/Box
PE02-10Fr	Foot Control 10Fr	20PCS/Box
PE02-12Fr	Foot Control 12Fr	20PCS/Box